#### Ionic Brushing Superpowers

# Ι N Se



Innovative cleaning through patented ionic technology.

Gentle-on-gums 31,000 strokes p/m sonic vibrations.

### 31,000 strokes p/m sonic vibrations

Our gentle-on-gums brushing has been specially designed to take better care of your mouth. There are 3 different brushing modes to suit your brushing needs.

## UV LED Light + patented Tio2 bar

UV light strikes TiO2 (Titanium Dioxide) bar and OH radical (hydroxyl radical) and IONs are generated.

Our superpowers lie in our patented technology.

Whereas most electric toothbrushes only remove some of the plaque and bacteria, ionic technology actually suppresses the bacteria causing plaque too.

# You can expect smoother, cleaner feeling teeth over time.

The OH radical and ION (OH-) generated by ION-Sei help to suppress bacteria attachment and therefore plaque formation.

# ION - Sei

# Weak current from NiMH (nickel-metal hydride) battery

The ION-Sei works with your body's natural electrical charge—producing an opposite weak charge that inhibits bacterial growth and helps to attract bacteria from teeth onto the brush.

Topi



#### ION-Sei specifications

#### Description

ION-Sei Sonic Ionic Electric Toothbrush

#### Rate voltage/Frequency AC100 - 240V 50/60Hz

AC100-240V 50/60HZ

#### **Battery type**

rechargeable NiMH (nickel-metal hydride)

#### **Charging time**

approx. 10 hours

#### Operating time

approx. 40 cycles (1 cycle = 2 min)

#### **Power consumption**

3W

#### Vibration speed

Sensitive 18,000 strokes/min Standard 24,000 strokes/min Deep 31,000 strokes/min

#### Size

diameter 31 mm, length 234 mm with brush, 210 mm without brush

#### Weight

96g\* with brush, 93g without brush \*Leading brand model approx. 135g

#### Waterproof

unit IPX5 charger stand IPX7

#### Parts included

ION-Sei handle Soft and standard brush USB cable EU and UK adaptor Charger stand Transparent charger cover.

#### **Brushing modes**

- Deep
- For a supercharged clean
- Standard
   For an everyday clean
- Sensitive For a gentle-on-gums clean



