Sonic electric toothbrush  
Elektrische Schallzahnbürste

# ION-Sei

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Read these brief operating instructions **before** first use of your sonic electric toothbrush and keep them for later reference.

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These operating instructions incorporate proprietary information which is subject to copyright law. No part of these operating instructions may be reproduced or transferred for any purposes whatsoever without prior written consent, independent of the means or of with which media this takes place: neither electronically nor mechanically.

Release 1, December 2017 – All rights reserved.

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## **These operating instructions are brief ...**

... and you only need a few minutes to read them. For the optimum use of your new ION-Sei sonic electric toothbrush, you need some information and must know a few simple steps. We will show you what and how on the following pages – as concisely as possible and in as much detail as required.

### **Intended use**

The ION-Sei sonic electric toothbrush is exclusively intended for cleaning and caring for human teeth. Any other use is expressly prohibited.

The ION-Sei sonic electric toothbrush is only designed for private use and not intended for application by several patients at a dentist's or a dental clinic. Its use is determined as contrary to the intended use if the device or parts thereof are not inserted as described in these operation instructions.

### **EU guidelines**

The ION-Sei sonic electric toothbrush complies with the requirements of the European guidelines 2014/30/EU (electromagnetic compatibility) and 2011/65/EU (restriction of the use of certain dangerous substances in electrical and electronic devices).

## Your safety is important!

Even if a sonic electric toothbrush is not dangerous, you should be aware of possible risks when handling it and act in a manner which neither endangers you nor other persons unnecessarily. Please therefore observe the following rules:

- » Only operate the sonic electric toothbrush in surroundings which are in accordance with the specifications in the section **Technical data** on page 14. Temperatures which are too high and excessive damp may damage the installed rechargeable battery and cause fire or an explosion. The supply voltage must be in accordance with the name plate specifications on the charger.
- » The charger is not waterproof. Do not use the charger outdoors. Do not touch the charger with wet hands. Ensure that the charger is always dry and cannot come into contact with water.
- » The sonic electric toothbrush is not a toy. Children from an age of 10 as well as persons with restricted physical, sensory or mental abilities or a lack of experience and knowledge may use the sonic electric toothbrush if under supervision or they have been instructed in its safe use and understand the dangers arising from it.
- » Do not use the sonic electric toothbrush if you have a pacemaker or other implants fitted.
- » If you are undergoing treatment for a dental disorder, consult your dentist before using the sonic electric toothbrush. If you experience excessive gum bleeding after use of the sonic electric toothbrush, seek advice from your dentist.
- » You cannot repair the sonic electric toothbrush, the charger and the USB cable yourself. If you find any damages, you must dispose of the damaged component in accordance with statutory regulations.

## Technology and function

In order to support optimum tooth cleaning and tooth hygiene, the sonic electric toothbrush ION-Sei combines high-frequency strokes with patented ion flow technology.

The sonic electric toothbrush head vibrates with a maximum high speed of 31,000 strokes per minute. You can choose between three vibration levels: sensitive, standard and deep.

The hydrodynamic effect arising from the vibration causes loosening, removal and de-struction of the bacterial biofilm on and between the teeth. As opposed to traditional toothbrushes, the brush head may only be pressed against the teeth lightly at an angle of 45 degrees due to the application of this technique. Due to the brush's own weight, the pressure applied is fully adequate. The low pressure makes dental hygiene considerably gentler than with a traditional toothbrush. Instead of making brushing movements, the brush head should be held as still as possible and moved from tooth to tooth to gain the best results.

If UV light (as a daylight component or from an installed UV lamp) is cast upon the moist titanium dioxide stick to which the brush head is attached, a flow of ions (electrically charged particles) is generated. The contact surface on the reverse of the sonic electric toothbrush ensures that the low-voltage current can flow over the surfaces of the teeth when the hand touches the contact surface, which is always the case if held normally by hand. The negatively charged ion flow electrons destroy the bacterial deposit and the bond between plaque and the teeth.


## Cleaning the teeth

Dentists recommend a cleaning time of at least two minutes. If using the sonic electronic toothbrush for the first time, it is possible that you experience a tickling or tingling sensation which is triggered by the high-speed vibrations. Following a brief period of use, this effect should no longer be felt.

You can use the sonic electric toothbrush with or without toothpaste. Special toothpaste is not required.


If you have not used a sonic electric toothbrush before, you should let your dentist or dental hygiene specialist show you how to use it correctly.

## Overview

- 1 Charging station with a transparent top
- 2 Handle
- 3 Brush head
- 4 USB cable
- 5 UK adapter
- 6 EU adapter

## Operating elements and displays

- 1 Titanium dioxide stick, for generating the ion flow
- 2 UV LED light (not visible), illuminated in blue when the toothbrush is switched on
- 3 Power switch for switching the toothbrush on and off and changing the vibration level
- 4 Three LEDs, display vibration level 1, 2 or 3 when the toothbrush is on
- 5 Charge control LEDs, display the current rechargeable battery's charge status (if full, three LEDs are illuminated)

# Preparations

## Unpacking

Open the packaging carefully and remove all the parts. Check all parts for complete and intact. If a part is visibly damaged, return everything complete to your dealer. Do not use any damaged parts.

Dispose of packaging not needed in accordance with statutory requirements.

## Setting up

Set up the charging station close to a mains socket on a stable and flat surface which it cannot fall off. Do not place the charging station on the edge of a tub or a sink. Ensure that the installation location is protected against splashing water and neither subjected to intense heat from radiators nor to strong sunlight.

Place the handle on the charging station and the brush heads on the two holders on the left and right behind the handle.

## Charging the rechargeable battery

Before initial use, the installed rechargeable battery must be fully charged. Full charging takes a maximum of 10 hours. Observe that the rechargeable battery is charged at normal room temperature. The ideal charging temperature is between 5 and 35 °C.

After initial full charging, you should charge the rechargeable battery roughly once per week or when the charge level indicates low.

- On the rear of the charging station, open the rubber cap above the connection socket.
- Pay attention to the correct alignment, plug the small USB cable plug carefully into this connection socket, and plug the large USB cable plug into the charger connection socket.
- Plug the charger into an appropriate mains socket.
- Place the handle on the charging station.

As soon as the charger is in the socket, the LED at the front of the charger illuminates. The three charge control LEDs at the front of the hand device illuminate alternately during the charging process. The charging process is completed when all three LEDs illuminate simultaneously.

- After charging, disconnect the charging device from the mains socket.
- Disconnect the small USB cable plug from the charging station connection socket and close the connection socket with the rubber cap in order to prevent the penetration of water.

For charging, you can use all devices with a Type A USB connection socket (USB 2.0 or higher) instead of the charger.

## Use

Dentists recommend being methodical when cleaning your teeth and to clean the four quadrants of the teeth (top right, top left, bottom left, bottom right) in sequence. In order to follow this procedure, the toothbrush is switched off briefly three times and then on again, respectively after 30 seconds. That is the signal to change quadrants. After two minutes (= 4 x 30 seconds), the toothbrush is completely switched off automatically. Naturally, you can switch it on again if you want to continue brushing.

To change the vibration level, press the power switch several times. Changes are always made in the same order: level 1 > level 2 > level 3 > off. Depending on the vibration level, one, two or three LEDs illuminate below the power switch. When the toothbrush is on, the blue UV LED light is also lit at the tip of the brush head.

**Attention!** Before attaching the brush head, pay attention to correct alignment. The bristles must point forwards towards the power switch on the handle. The little ridge at the top of the handle must be guided into the groove at the bottom of the brush head.

- Select a (soft or hard) brush, align it, guide it over the titanium dioxide stick at the top of the handle and attach it firmly as far as it will go.
- Dampen the bristles under flowing water and apply toothpaste to the bristles as desired.
- Hold the handle in one hand and press the brush head lightly against your teeth, best at an angle of 45 degrees to the gum.
- Press the power switch once or several times briefly using your thumb to switch on the toothbrush and set the desired vibration level.
- If you press the power switch again briefly after activating the third level or keep it pressed longer, the toothbrush is switched off again.
- Guide the bristles from one tooth to the next, both from the front and back with a small back and forth motion.
- During the brief automatic interruption, change the quadrant.

When the toothbrush has switched off automatically after two minutes, it starts at the previously selected vibration level when switched on again.

## Cleaning and maintenance

**Attention!** For cleaning, do not use detergents, alcohol, solvents and also no abrasive or hard items. They can seriously damage the surfaces. Do not clean any parts in the dish-washer.

- » Rinse the brush head thoroughly under flowing water after cleaning your teeth and then shake out the water.
- » Remove the brush head at least once per week and clean the contact surfaces and titanium dioxide stick with a soft cloth.
- » Clean all the flat surfaces regularly with a dampened, soft cloth.
- » Pull the transparent charging station top off upwards to clean both parts of the charging station separately. Place the transparent top back on when ready.
- » Do not use brush heads with bent or crushed bristles. Replace the brush head at the latest after three months or as soon as they show signs of wear.
- » Only use the brush heads recommended by the manufacturer for this sonic electric toothbrush.
- » Place unused brush heads on the charging station holder so that they can dry and the bristles do not come into contact with impurities.
- » Ensure that no body care products (for example hand lotion, hair gel or oil) come into contact with the sonic electric toothbrush, the charger or the cable. These products can damage the surface and lead to cracks or breaks.

## Passing on / disposal

Only pass on the ION-Sei sonic electric toothbrush together with these operating instructions.

The sonic electric toothbrush and the charger contain high-quality materials and components which can be recycled and reused.

The symbol means that an electric or electronic device must be disposed of separately from normal household waste. In the European Union, there are various collection systems for used electric and electronic products. Dispose of the devices marked like this at your local collection and recycling point.

Before disposing of the handle, you must remove the rechargeable battery from it and dispose of it separately. The rechargeable battery must not be disposed of with household waste.

If you are unsure or do not have the tools required, have the rechargeable battery removed by your dealer or a specialist.

**Attention!** This process is irreversible. If the cap at the base of the handle has been removed, the toothbrush may not be used.

- Remove the toothbrush from the charger cradle, switch it on and let it run until the rechargeable battery is completely empty.
- Unscrew the screw in the base of the handle with a Philips type screwdriver and remove the cap.
- Loosen the glued on spiral wire and bend it to the side.
- Push a small screwdriver into the gaps and release the two locks between the interior part and the casing by rotating it slightly.
- Hold the handle upside down and press the titanium dioxide stick upright on a table top to press the inner part out of the casing on the open side.
- Remove the rechargeable battery and the cable plug.

## Troubleshooting and technical support

The following table is intended to help you to identify the cause of any malfunctions and, as applicable, to resolve them.

If you are unable to rectify a malfunction, please refer to your specialist dealer or our customer hotline (see page 3).

<b>Malfunction</b>	<b>Possible cause</b>	<b>Recommended remedy</b>
The toothbrush cannot be switched on.	Rechargeable battery empty.	Charge the rechargeable battery.
Unaccustomed loud noises after switching on.	Brush head not attached properly.	Remove brush head, align it correctly and attach firmly as far as it will go.
Tingling sensation during use.	Normal reaction during initial use.	Following a brief period of use, this effect should no longer be felt.
Short interruptions whilst cleaning your teeth.	A normal automatic function, information to change quadrants (see page 10).	
Rechargeable battery is not charged.	No charging current.	Check whether both USB cable plugs are mechanically secured.  Plug the charger into the socket.  If the USB socket from a different device is used, switch the device on.
Charger LED is illuminated but the rechargeable battery is not charged.	Charger or rechargeable battery faulty.	Have the charger checked by a specialist dealer.

## Technical data

Nominal charger voltage	100 – 240 V ~, 50/60 Hz, 0.2 A
Rechargeable battery type	3.7 V NiMH (nickel-metal hydride)
Rechargeable battery life	Approx. 40 cycles, each 2 minutes
Recharging time	Max. 10 hours
Power consumption	3 W
Vibration speed	3 levels: Sensitive (18,000 strokes/min.) Standard (24,000 strokes/min.) Deep (31.000 strokes/min.)
Connections	Charging station: Micro USB socket for connection of the USB cable for charging the rechargeable battery Charger: USB type A socket for connection of the USB cable to the charging station
Environmental conditions	Operation: 5 – 35 °C, 80 % rel. humidity Storage: -10 – 40 °C, 80 % rel. humidity
Handle dimensions	Length: 210 mm without brush head 234 mm with brush head  Diameter: approx. 31 mm
Handle weight	96 g with brush head
Protection class	Handle: IPX5 Charging station: IPX7 (USB connection socket closed)

## Further information options

If you have any questions regarding the use of the ION-Sei sonic electric toothbrush which are not covered by these operating instructions, please refer to your dealer.

Further current information regarding the ION-Sei sonic electric toothbrush is available on the Internet on the **[www.ion-sei.com](http://www.ion-sei.com)** pages.

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